# BEAUTIFLY

Scrub Mini

Young  
series

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## Product parameters:

<b>Rated Voltage:</b>	AC 100V~240V 50HZ/60HZ
<b>Output Voltage:</b>	5V/400mA
<b>Battery:</b>	Lithium Battery
<b>Net Weight:</b>	70g
<b>Ultrasounds:</b>	25KHZ
<b>Power Supply:</b>	2W
<b>Unit Size:</b>	140*42*20mm
<b>Set contains:</b>	manual, USB cable, other accessories

## Safety statement

Please read the instructions carefully before use. Follow these safety guidelines when using device.

## Part description

1. Stainless steel spatula
2. Protection cap
3. Power button/Mode key
4. USB charging interface

## Product Functions

- Ensures deep cleaning of the skin.
- Reduces color spots and decomposes pigment.
- Minimizes deep wrinkles and improves skin elasticity.
- Promotes blood circulation and lymphoid metabolism.
- Enhances the regeneration ability and permeability of skin cells.
- Improves condition and appearance of the skin
- Is easy to carry and navigate.

## Charging

Fully charge the device before its first use.

- Check the status of the lithium battery before charging.
- Connect the USB charging cable to the USB interface.
- Connect to the USB interface of your computer, adapter (sold separately) or lithium battery of your cell phone.
- The LED flashes continuously during charging, and lights up when charging is complete.

## Charging Attentions

- Charge the device if it does not display the boot status when long press the power button.
- In operation, the device will make a beep-beep warning sound, then the shutdown state will appear. It indicates that the battery's margin is insufficient. Charge the product immediately.
- After working for 10 minutes, the product automatically goes into standby mode, and if not used for one minute, it turns off automatically.
- It takes approximately 120 minutes to fully charge. (Before first use, charging takes more than 5 hours.)
- The device can be used for about 60 minutes when the battery is fully charged.
- It is recommended to charge the device periodically.

## Operation guide

The switch lights flash when the mode is changed.

### Cleansing mode

1. Long press the power button for 1.5 seconds to turn on the device.
2. The first Cleansing mode will start automatically after the device is powered on. Daily care: reduction of blackheads, whiteheads, and acne.
3. Ultrasonic vibrations penetrate deep into the muscles, deeply cleansing the skin. Use with the facial cleanser.

## Moisturizing mode

1. Long press the power button for 1.5 seconds to turn on the device.
2. Press the power button one time to start the second mode, which is Moisturizing.
3. The ultrasonic vibrations and the negative ion make the active substances in cosmetics penetrate deep into the skin.

## Lifting Mode

1. Long press the power button for 1.5 seconds to turn on the machine
2. Short press the power button 2 times to start the third mode - Lifting.
3. Microcurrents generated by ultrasound and EMS stimulate cell activity. The process stimulates the production of collagen and elastin, which makes the skin toned and smooth.

## Cleaning and maintenance

- Protect the instrument from moisture, and do not place it in a dishwasher or microwave.
- Turn off the machine before cleaning.
- Do not use sharp, corrosive, or volatile thinner cleaning materials.
- Use a soft cloth or cotton pad moistened with cold or warm water to wipe gently.
- Do not rinse the device directly with water.

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Electrical and electronic devices as well as batteries and accumulators must not be mixed with other household waste. This is indicated by the symbol of the crossed out basket on the packaging. Worn-out electric and electronic equipment as well as batteries and accumulators should be delivered to local separate waste collection points or to the seller, and for details inquire in your commune. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage to them may lead to health detriment or environmental pollution during further disposal or recycling of the equipment.

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The importer hereby declares, at their sole responsibility, that the Scrub Mini Young product complies with the Directive of the European Parliament and of the Council 2014/30/EU of 26 February 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, as well as the Directive of the European Parliament and of the Council 2011/65/EU of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU of the European Parliament and of the Council as regards the list of restricted substances. The full text of the declaration of conformity is available at the website address.

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